



cryo:one -85°C

The Evolution in Whole Body Cryotherapy

The **cryo:one** from the market leader MECOTEC is the best and smallest full electric cold chamber in the world. Its compact and highly efficient room size enables you to experience all the advantages of the whole body cryotherapy (WBC) down to -85°C even in the smallest of spaces.

Instead of using liquid nitrogen, all of MECOTEC's cryo chambers cool with electric compressors. This allows you to safely enjoy your whole body cryotherapy in fresh air and with full vision and freedom of movement. Your entire body benefits 100% from the cold – from head to toe.



Redefine yourself with
cryo:one -85°C

With a duration of just 3 minutes,
whole body cryotherapy makes
long-lasting benefits possible:

- Faster recovery after training
- Enhanced rehabilitation after (sporting) injuries
- Rejuvenation of body, skin and mind
- Stimulating of metabolism



Starting from private use via usage in practises, right to commercial use in studios and spas – the **cryo:one** adapts to your needs. Customize and colorize the cold chamber according to your own ideas.



Benefits of whole body cryotherapy

Cold is your efficient and healthy way to physical and mental balance. At -85°C, the air contains no moisture, that's why the extreme cold feels much more pleasant. Due to the volume of air shrunk coming from these temperatures, you take in a much higher amount of oxygen with each breath. During the application, the body releases endorphins and dopamine and you feel energized immediately after the application – with simultaneous relaxation. The whole body cryotherapy in MECOTEC's electric chambers offers you excellent benefits!

- **Pain relief**
- **Faster regeneration**
- **Physical and mental well-being**
- **Strengthening the immune system**



Professional Refrigeration Technology – Made in Germany

The installation of the **cryo:one** is possible in rooms with a ceiling height from 2.50 m. For commissioning, a power and water connection are required. Our professional service technicians ensure that your **cryo:one** is expertly commissioned and ready for use within a few hours. The cold chamber can also be easily controlled via smartphone app or tablet.

We are the founders of commercial whole body cryotherapy with over 20 years of experience in developing our products.

MECOTEC stands for technology, production, installation and service from a single source in the field of cryotherapy. With currently more than 450 cryotherapy chambers worldwide.

Invest in the original.

Technical data

Chamber outside (Length x Width x Height)	1,2 x 1,2 x 2,5 m
Chamber inside (Length x Width x Height)	0,9 x 0,9 x 2,1 m
Engine (Length x Width x Height)	1,0 x 0,6 x 1,6 m
Weight of the whole unit	590 kg
Connection	400V AC +/- 10%, 3 Ph/N/Earth, 50/60 Hz, 25 A
Condenser	water-cooled (air-cooling on request)

MECOTEC
FREEZING TECHNOLOGY

We are your partner for cryotherapy.

www.mecotec.net

www.cryoone.net